

Bending expectations of conformability and stability.

Enhanced conformability


The hybrid cell structure is designed to enhance stent opening and conformability in bifurcations and tight curves.

Ease of use

All sizes of the Neuroform Atlas Stent are deliverable through Excelsior® SL-10® and Excelsior XT-17™ Microcatheters.

Higher deployment accuracy

The Neuroform Atlas Stent is designed to have very low foreshortening, which enables very high deployment accuracy.



Neuroform Atlas® STENT SYSTEM

The Neuroform Atlas Stent System is authorized under a Humanitarian Device Exemption (HDE). IRB approval is required prior to use.

Copyright © 2019 Stryker
AP002413 v1.0 | Page 1 of 2